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
Taq DNA Polymerase

The Polymerase Chain Reaction (PCR) is a powerful and sensitive technique for DNA amplification. Taq DNA Polymerase is an enzyme widely used in PCR. The following guidelines are provided to ensure successful PCR using NEB's Taq DNA Polymerase. These guidelines cover routine PCR. Amplification of templates with high GC content, high secondary structure, low template concentrations, or amplicons greater than 5 kb may require further optimization.

Step	Temperature	Time
1. Initial Denaturation at 95°C	95°C	30 sec
2. 30 cycles of: Denaturation at 95°C	95°C	30 sec
3. 30 cycles of: Annealing at 55°C	55°C	30 sec
4. 30 cycles of: Extension at 72°C	72°C	30 sec

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

The screenshot displays the Labguru dashboard with the following components:

- Header:** Includes the Inserm clé logo, the Labguru brand name, a support dropdown, and a user profile icon labeled "DSI Lab".
- Navigation Bar:** Features tabs for "e-Notebook", "Knowledgebase", "Inventory", and "Storage & Equipment". A "Recently Viewed" section with a search bar is also present.
- Alerts Section:** A row of cards titled "Actions require your attention" (58 total). The cards include:
 - "Signed experiments need witnessing" (0) with a "Go to experiments" button.
 - "Low stock alerts" (0) with a "Go to shopping list" button.
 - "Overdue tasks" (5) with a "Go to tasks" button.
 - "Unattached images & attachments" (53) with "Go to images" and "Go to attachments" buttons.
 - "Expired Stocks Alert" (54) with a "View Expired Stocks" button.
 - "Equipment require calibration" (1) with a "Go to equipment" button.
- Navigation Tabs:** "Dashboard" (active), "Activities", "Calendar", "Whiteboard", and "Latest Results".
- Calendar View:** Shows the date "October 26, 2018". The main content area displays a "Welcome to Labguru!" message and a list of recent activities:
 - Today at 15:23: Maria Lerner - Added a new experiment ELISA to folder New Folder in project My IPO Project.
 - Today at 15:23: Daria Inbar - Added a new experiment IPO Activation to folder My Experiments in project - 2018-03-04.
 - Tuesday at 15:21: Daria Inbar - Updated the bioinformatics.
 - Sunday at 11:02: Daria Inbar - Added a new project My project - 2018-03-04.
- Calendar Widget:** A weekly view for Friday, showing a time slot from 11:00 to 12:00 labeled "Séminaire".
- Annotations:** A purple diagonal watermark "Exemple-" is overlaid on the left side. A purple speech bubble in the bottom right corner says "Accès au support" with a support icon.

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Support ▾
DSI Lab ▾

e-Notebook
Knowledgebase
Inventory
Storage & Equipment
Recently Viewed

Mass spec for media protein
Mass spec for media protein
Mass spec pilot

Pc and L1 30 sample preparation and calibration (duplicate)

🔗 This experiment is a duplicate of Pc and L1 30 sample preparation and calibration

👤 Paul-Guy Dupré
👤 You and 37 more can edit
🚩 Flag Experiment
👤 Sign
More ▾

📅 Start date: Jan 19, 2018
⌚ Duration (days): Enter duration

🏷 tag this item

Description

🔗
💬
More ▾
📄 Copy to report
Save

Instructions are provided in this section to use the NeonTM device with the NeonTM Pipette Station and NeonTM Kits for electroporation of mammalian cells.


General instructions to prepare cells for use with the NeonTM Transfection System are described below. For primary and stem cell types, use the established methods developed in the laboratory.

Table 1

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Inserm clé labguru Support ▾ DSI Lab ▾

e-Notebook Knowledgebase Inventory Storage & Equipment Recently Viewed Search...

Protocols

PureYield™ Plasmid Maxiprep System

Promega You and 37 more can edit

⌚ Duration (days): Enter duration

🏷 tag this item

Description

As research moves from DNA sequencing to analysis of gene function, the need has increased for rapid methods by which to isolate large quantities of high-quality plasmid DNA. The PureYield™ Plasmid Maxiprep System is designed to isolate high-quality plasmid DNA from eukaryotic and prokaryotic cells. The system provides a rapid method for purification using a silica-membrane column. Plasmid DNA can be purified in approximately 60 minutes, greatly reducing the time spent on purification compared to other membrane column methods. The PureYield™ Plasmid Maxiprep System also incorporates a unique Endotoxin Removal Wash designed to remove substantial amounts of protein, RNA and endotoxin contaminants from purified plasmid DNA, and is suitable for sensitive applications such as eukaryotic transfection, in vitro transcription and cell-free expression. Purification is achieved without isopropanol precipitation of purified plasmid DNA, providing rapid purification as well as a high concentration of pure plasmid DNA. The PureYield™ Plasmid Maxiprep System purifies up to 1mg of plasmid DNA sensitive applications such as eukaryotic transfection, in vitro transcription and cell-free expression. Purification is achieved without isopropanol precipitation of purified plasmid DNA, providing rapid purification as well as a high concentration of pure plasmid DNA. The PureYield™ Plasmid Maxiprep System purifies up to 1mg of plasmid DNA with an $A_{260}/A_{280} > 1.7$ from 250ml of overnight bacterial culture, transformed with a high-copy-number plasmid. The system requires a vacuum pump and manifold (e.g., the Vac-Man® Laboratory Vacuum Manifold, 20-sample [Cat.# A7231]), a centrifuge with a fixed-angle rotor for lysate clearing and a tabletop centrifuge with a swinging bucket rotor for elution by centrifugation for optimal DNA recovery and yield, elution can be performed by vacuum, with the Eluator™ Vacuum Elution Device (Cat.# A1071; see Section 4.B). Use of the Eluator™ Elution Device eliminates the need for a tabletop centrifuge and swinging bucket rotor.

For Product Components and Storage Conditions please [click here](#)

add a section

Preparation and Lysis of Bacterial Cell Cultures

Note:

Throughout the remainder of this document the supplied Cell Resuspension Solution (CRA), Cell Lysis Solution (CLA) and Neutralization Solution (NSB) are referred to as Cell Resuspension Solution, Cell Lysis Solution and Neutralization Solution, respectively

Note:

Perform all lysis steps at room temperature

Steps


1	Grow 100-250ml of transformed E. coli bacterial cell culture overnight (16-21 hours) at optimal culture conditions.	00 : 00 : 00 Completed SMS
2	Pellet the cells by centrifugation at 5,000 xg for 10 minutes at room temperature, and discard supernatant. Drain tubes on a paper towel to remove excess media	00 : 10 : 00 Completed SMS

☐ All steps completed

Linked resources

- Duplicate...
- Archive
- Sign
- Linked resources
- Editable by...
- Visible to...
- Share to repository
- Version history
- Print
- Delete protocol...

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
The screenshot displays the 'labguru' inventory management system interface. The top navigation bar includes the 'Inserm clé' logo, the 'labguru' name, and a 'Support' dropdown. Below this, a secondary navigation bar lists 'e-Notebook', 'Knowledgebase', 'Inventory', and 'Storage & Equipment'. The main content area features a header with '58 Actions require your attention' and a 'See all' link. Below the header, there are four alert boxes:

- Signed experiments need witnessing** (0 alerts): Includes a 'Go to experiments' button.
- Expired Stocks Alert** (54 alerts): Includes a 'View Expired Stocks' button.
- Low stock alerts** (0 alerts): Includes a 'Go to shopping list' button.
- Equipment require calibration** (1 alert): Includes a 'Go to equipment' button.

Below the alerts, there is a grid of equipment cards. Each card displays a number in parentheses (e.g., 6 (A6), 7 (A7), etc.) and a green checkmark icon. Three cards in the middle row are highlighted with a blue border and contain the text 'Nuclease ... 1/12/18'.

Evolution des instruments: la calibration

The screenshot displays the 'labguru' web application interface. The top navigation bar includes the 'Inserm clé' logo, the 'labguru' brand name, and a 'Support' dropdown menu. Below this, a secondary navigation bar lists 'e-Notebook', 'Knowledgebase', 'Inventory', and 'Storage & Equipment'. The main content area shows a breadcrumb trail 'Equipment > centrifugeuse 2' and the title 'centrifugeuse 2'. A tabbed interface at the bottom of the main content area includes 'Info', 'Experiments & Protocols', 'Output Files', and 'Calibrations'. The 'Calibrations' tab is active, displaying 'centrifugeuse 2 Calibrations'. On the right side of this tab, there is a blue 'Calibrate' button with a scale icon. Below the title, a table with three columns is visible: 'Calibration data', 'Performed by', and 'Performed at'. The 'Performed by' column contains the name 'Sylvain Demey' and the 'Performed at' column contains the timestamp 'Oct 05, 2018 11:03'. The 'Calibration data' column is currently empty.

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